

Lyncmed Product Profile



Global Linker, Healthcare Better; Dream Maker, Win Together >>>



Disposable Face Mask Product Profile

KA International Business Dept

Lyncmed Medical Group

Face Mask Packing & Artwork



New Box Design Artwork



Previous LM Face Mask Box Picture



Face Mask Carton Design Artwork



Face Mask Carton Design Artwork

Face Mask Certificates –EN14683 BEF Test Report



Nelson
 Medical Technology (Beijing) Co., Ltd.
 Room 115, No. 1111
 South Huata Road, Chongqing District
 Beijing 100020
 China

Bacterial Filtration Efficiency (BEF) and Differential Pressure (Delta P) Final Report

Test Article: Medical Mask Non-woven Face mask
Lot No.: 20200716
Manufacturer: 23 Aug 2019
Study Number: 100011001
Test Facility: Nelson Laboratories, LLC
 8280 S. Redwood Rd.
 San Jose, CA 95128 U.S.A.
Test Procedure: Standard Test Protocol (STP) Number: STP0008 Rev. 10
Traceability: None

Summary: The BEF test is performed to determine the filtration efficiency of test articles for comparing the technical control counts upstream of the test article to the technical control downstream. A suspension of *Staphylococcus aureus* was aerosolized using a nebulizer and delivered to the test article at a constant flow rate and tested as described. The challenge volume was measured as 17.2 ± 0.1 mL using a gravimetric method with a mean particle size (MPPD) of 2.2 ± 0.3 µm. The particles were drawn through a 90 mm, 0.45 µm pore size, membrane filter for collection. The test method conforms with EN14683:2019, Annex B, and Annex D(2)(3).

The Delta P test is performed to determine the breathability of test articles for measuring the differential air pressure on either side of the test article using a manometer, at a constant flow rate. The Delta P test was designed to comply with ISO 9001:2015, Section 4.4.1 and conforms with EN 14683:2019, Annex C, and Annex D(2)(3).

All test method acceptance criteria were met. Testing was performed in accordance with US FDA good manufacturing practice (GMP) regulations 21 CFR Parts 312, 314, and 315.

Test Data: None
BEF Test Area: 48 cm²
BEF Flow Rate: 20 L/min per minute (2.0/min)
Delta P Flow Rate: 3 L/min
Compliance Requirements: BE ≥ 95% media permeability (95% and 21 ± 2%) for a minimum of 4 hours
Test Article Dimensions: 177 mm x 108 mm
Positive Control Average: 0.0 ± 0% CPV
Negative Monitor Count: 0 CPV
MPD: 2.2 µm



Issued by: *Janella Pong* Issued on: 2020-08-02 Issued for: *Medical*
 Nelson Laboratories, LLC 8280 S. Redwood Rd. San Jose, CA 95128 U.S.A.



Nelson
 Medical Technology (Beijing) Co., Ltd.
 Room 115, No. 1111
 South Huata Road, Chongqing District
 Beijing 100020
 China

Results

Test Article	BEF (%)	Delta P (mmHg)	Media Permeability (%)
1	99.0	1.0	92.1
2	99.0	1.0	92.0
3	99.7	1.0	92.0
4	99.0*	1.0	92.0
5	99.0	1.0	92.0

* Values were not detected either on any of the Anderson sampler cups for the test article

The BEF efficiency percentage were calculated using the following equation:

$$BEF = \frac{C - F}{C} \times 100$$

C = Positive control average
 F = Filter count (total measured concentration of the test article)
 Note: The filter count total is available upon request

Lyncmed Company Certificate



ISO-13485 Certificate

3 Layer Non-woven Disposable Face Mask

Disposable Face Mask Key Features:

- Skin Friendly High Quality PP Material, 3-Ply
- Low Breathing Resistance, Bacterial Filtration Efficiency(BEF)>98%
- Ear Loop, Elastic Band, Latex Free
- Anatomic Adjustable Integrated nose bridge
- Size:17.5*9.5cm



Face Mask CE Certificate

EC DECLARATION OF CONFORMITY

Manufacturer: Liyemed Medical Technology(Beijing) Co., Ltd.
Room 110, No.1111 South Hualu Rd, Changyang District 10002 Beijing, China

whose single Authorized Representative: BSI(UK) Certification Company Limited
8001, MAPLE HOUSE, 118 HIGH STREET, PURLEY, LONDON, ENGLAND

We declare under our sole responsibility that: Also covers products: Surgical gowns, Surgical drapes, Surgical packs, Prosthetic Caps, Coveralls, Face masks, Goggles, Underwear(Gloves), Goggles, Goggles, Masks, Pads, Shoe covers, Aprons, Gowns, Trousers, Sterile covers, Floor covers, Beard covers, Bedgowns, Suits, Suits, Hair Suits, Disposable PE products: Shoe covers, Tey covers, Footwear, Caps, Goggles, Aprons, Gowns, Shoe covers, Slips, Sheets, Pillow Case, Pad covers, Sides, Disposable Vinyl examination gloves, Disposable Nitrile examination gloves, Disposable Latex examination gloves, Latex Glove for single use, Latex Glove for single use

the medical device: None

of class: Class I

under the provisions of the Directive 93/42/EEC (or 90/269/EEC) which apply to it: Medical Device Directive 93/42/EEC

Conformity assessment procedure: EN45001, Annex Y

The above mentioned declaration of conformity is necessary under the responsibility of: Liyemed Medical Technology(Beijing) Co., Ltd.

Place, date: Beijing, China, Management Representative
Beijing 09th Dec. 2017 Legally binding signature

(Red circular stamp with a star and Chinese characters)

(Handwritten signature)

EC Declaration of Conformity

CE Certificate

Face Mask Certificates –EN14683 SBPR & EN14683 Blood Penetration Resistance Report



Nelson
 Lynwood Medical Technical (Beijing) Co., Ltd
 Room 111, No. 1111,
 South Huhai Road, Chiyang District,
 Beijing 100020,
 CHINA

Synthetic Blood Penetration Resistance Final Report

Test Article: ProJet Normal Non-woven Face Mask
 LPT #20244719
 Study Number: 1089212-001
 Study Planned Date: 23 Aug 2018
 Testing Facility: Nelson Laboratories LLC
 4380 S. Redwood Rd
 Red Lake City, UT 84003 U.S.A.
 Test Procedure(s): Standard Test Protocol (STP) Number: STP0012 Rev 08
 Deviation(s): None

Summary: This procedure was performed to evaluate surgical face masks and other types of protective clothing materials designed to protect against fluid penetration. The purpose of this procedure is to simulate an adverse spray and evaluate the effectiveness of the test article in protecting the user from possible exposure to blood and other body fluids. The distance from the target area surface to the tip of the particle is 30.5 cm. A test volume of 2 mL of synthetic blood was employed using the targeting area method.

This test method was designed to comply with ASTM F1982 and ISO 22820 (as referenced in EN 14683 (2016) and GB26263 (2016)) with the following exception: ISO 22820 requires testing to be performed in an environment with a temperature of 21 ± 0.5°C and a relative humidity of 50 ± 10%. Instead, testing was performed in ambient conditions within one minute of removal from the environmental chamber used at these parameters.

All test method acceptance criteria were met. Testing was performed in compliance with US FDA good manufacturing practice (GMP) regulations 21 CFR Parts 210, 211 and 820.

Number of Test Articles Tested: 32
 Number of Test Articles Passed: 32
 Test Site: Outside
 Pre-Conditioning: Minimum of 4 hours at 21 ± 0.5°C and 50 ± 10% relative humidity (RH)
 Test Conditions: 18 ± 0.5°C and 50% RH

Results: Per ASTM F1982 and ISO 22820, an acceptable quality level of 0% is met for a normal single spraying plan when 420 of 32 test articles show passing results.

Test Pressure: 100 mmHg (10.0 kPa)	Particle Blood Penetration
1-18, 20-26, 28-32	None Seen
19, 27	Yes

Study Director: *Janeille King*
 Study Engineer: *Brandon L. Williams*
 Study Completion Date: *10/2/2018*



Sponsor
 Dinnia Dong
 Lynwood Group
 No. 1111, South Huhai Road,
 100020, Beijing CHINA

EN 14683 2005 Synthetic Blood Penetration Resistance Final Report

Test Article: Non-woven Face Mask
 Lab Article Number: F00304
 Study Received Date: 27 Jan 2018
 Test Procedure(s): Standard Test Protocol (STP) Number: STP0012 Rev-05

Summary: This procedure was performed to evaluate surgical face masks and other types of protective clothing materials designed to protect against fluid penetration. The purpose of this procedure is to simulate an adverse spray and evaluate the effectiveness of the material in protecting the user from possible exposure to blood and other body fluids. The distance from the target area surface to the tip of the particle is 30.5 cm. A test volume of 2 mL of synthetic blood was employed. The test method was designed to comply with ASTM F1982 and EN 14683 2005. All test method acceptance criteria were met.

Number of Test Articles Tested: 32
 Number of Test Articles Passed: 32
 Test Site: Outside
 Pre-Conditioning: Minimum of 4 hours at 21 ± 0.5°C and 50 ± 10% relative humidity (RH)
 Test Conditions: 22 ± 0.5°C and 20% RH

Study Director: *Brandon L. Williams*
 Study Engineer: *Brandon L. Williams*
 Study Completion Date: *10/2/2018*



Face Mask Certificates –EN14683 MC Test Report



Nelson Labs.
A Sirona Health company

Sponsor:
Masks Co.
Lynwood Medical Technical (Shang) Co.
Room 118, No. 1111, South Huifu Rd., Changyang District,
Beijing 102200,
China

Microbial Cleanliness (Bioburden) of Medical Masks Final Report

Test Article: Product Name: Non-woven Face mask
SKU: MC044714
Study Number: 138874-001
Study Protocol Date: 22 Aug 2018
Testing Facility: Nelson Laboratories, LLC
6281 S. Redwood Ave.
Salt Lake City, UT 84123 U.S.A.
Test Procedure(s): Standard Test Protocol (STP) Number: STP00308 Rev. 14
Customer Specification Sheet (CSS) Number: 201900308 Rev. 01
Deviation(s): None

Summary: This testing was conducted in accordance with EN 14683:2014, with the exception of appropriate volumes of eluent used when performing the extraction procedure and a temperature range of 30-35°C used for aerobic incubation.
When bioburden results are calculated using a software program, manual calculations may differ slightly due to rounding. The counts determined on products are colony forming units and may not always reflect individual microorganisms. The sponsor performs an optical analysis and determines the acceptable limits. Testing was performed in compliance with US FDA good manufacturing practice (GMP) regulations 21 CFR Parts 201, 201.1 and 625.

Unit Number	Weight (g)	Sample	T/Fungal	Total Aerobic CFU/mL	Total Bioburden CFU/mL
1	3.2	88	12*	88.7	21.1
2	3.3	48	2*	51.8	16.5
3	3.3	33	3*	35.0	10.9
4	3.4	51	<5	54.2	16.9
5	3.5	99	2*	77.8	22.8

Necessary Efficiency: 85.7%

* = No Organisms Detected

Note: The results are reported as colony forming units (CFU) per mask.
Note: Sample positive testing was performed using *Bacillus anthracis*. The test article was not positive using this test method.

* Spread: Count is considered a minimum estimate due to swarming of certain colonies on the membrane.

Study Director: 
Nelson Laboratories, LLC
Study Completion Date:  

480.986.7666 | 480.986.7666 | 480.986.7666 | 480.986.7666
Page 1 of 2



Nelson Labs. Microbial Cleanliness (Bioburden) of Medical Masks Final Report

Study Number 138874-001

Test Method Acceptance Criteria: If applicable, aseptic controls are acceptable for the bioburden test results. The number of masks to be tested shall be a minimum of 5 or more to meet an acceptable quality level of 0%. The bioburden of the medical mask shall be <= 25 CFU/g mask.

Procedure:

Positive Controls/Membranes:	<i>Bacillus anthracis</i>
Control Fluid:	Phosphate Tryptone (PT) Solution Chloride
Extract Fluid Volume:	~200 mL
Extract Method:	Control Swabbing for 5 minutes at 230 rpm
Filtering Method:	Membrane Filtration
Agar Medium:	Tryptic Soy Agar
Incubation (Bacteria):	Sabouraud Dextrose Agar with Chloramphenicol
Extraction Time Method:	Extraction Time Method
Aerobic Bacteria:	Plates were incubated 3 days at 30-35°C, then enumerated
Fungal:	Plates were incubated 7 days at 30-35°C, then enumerated

480.986.7666 | 480.986.7666 | 480.986.7666 | 480.986.7666
Page 2 of 2